




FreeStyle  
Libre 14 day  
FLASH GLUCOSE MONITORING SYSTEM

## Sharing Your Reports

LibreView allows you to share glucose data with your healthcare providers

If you are uploading your reader data at home or using FreeStyle LibreLink app, your provider may ask you to share your glucose data with them.



Download the FreeStyle LibreLink app on your iPhone\* from Apple App Store. The FreeStyle LibreLink app is compatible with iPhone 7 and later running iOS 11 and later.

There are two ways to connect your account to your physician. First, create your LibreView account at [www.LibreView.com](http://www.LibreView.com) or through FreeStyle LibreLink app, then follow the directions below:

### Option 1: Accept the invitation from your healthcare provider

- Your provider will ask for a current email address
- You will receive an email invitation to that email address
- Open the email and click “Accept”
- Sign in to your LibreView account and click “Accept”

### Option 2: Log in manually with your healthcare provider’s ID

- Your doctor will provide you with a “Practice ID” 122M4834
- Log in to [www.LibreView.com](http://www.LibreView.com) and navigate to “Account Settings” from the right navigation menu
- From the “Account Settings” page, click on the “My Practices” tab
- Input the “Practice ID” and click “Add”

**Need Help?** Contact LibreView support at 855-786-4263 or visit the “Support” section on [www.LibreView.com](http://www.LibreView.com)



life. to the fullest.™

Abbott

\*The FreeStyle LibreLink app is compatible with iPhone 7 and later running iOS 11 and later.

LibreView is developed, distributed, and supported by NewYu, Inc. The LibreView data management software is intended for use by both patients and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis and evaluation of historical glucose meter data to support effective diabetes management. The LibreView software is not intended to provide treatment decisions or to be used as a substitute for professional healthcare advice. Trademarks are the property of their respective owners. Images are for illustration purposes only. Not actual patient data.